## Lavori scientifici / Clinical evidences

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### **The kalos device: experimental protocol for airway clearance in patients with ineffective cough**

S. Belli<sup>1</sup>, F. Masocco<sup>1</sup>, D. Cattaneo<sup>1</sup>, G. Savio<sup>1</sup>, B. Balbi<sup>1</sup>  
<sup>1</sup> ICS Maugeri, Istituto Scientifico di Veruno (Veruno, IT)  
5<sup>th</sup> ERCA Lione 2018

### **A pilot study on the non-invasive management of tracheobronchial secretions in tracheotomised patients**

S. Belli<sup>1</sup>, D. Cattaneo<sup>1</sup>, F. D'Abrosca<sup>2</sup>, I. Prince<sup>1</sup>, G. Savio<sup>1</sup>, B. Balbi<sup>1</sup>.  
<sup>1</sup> ICS Maugeri, Istituto Scientifico di Veruno (Veruno, IT); <sup>2</sup> IRCCS Ca' Granda-Ospedale Maggiore Policlinico (Milano, IT)  
5<sup>th</sup> ERCA Lione 2018

### **Management of the airway in the tracheostomy patient in critical areas: a comparison of two devices for a pilot study.**

L. Cavanna<sup>1</sup> PT, C. Bardonè<sup>2</sup> PT, M. Bertoni MD<sup>3</sup>, P. Severgnini MD<sup>4</sup>  
<sup>1,3</sup> Rehabilitation Department ASST-Sette Laghi, Varese; <sup>2</sup> Pulmonology Department ASST-Sette Laghi, Varese; <sup>4</sup> University of Insubria - DBSV – ICU ASST Settelaghi Varese, Italy;  
5<sup>th</sup> ERCA Lione 2018

### **Running title: EFA and pep physiotherapy after pulmonary lobectomy**

N. Rotolo<sup>1</sup>, M. Cattoni<sup>1</sup>, M. D'Andria<sup>1</sup>, L. Cavanna<sup>1</sup>, G. Patrizio<sup>1</sup>, C. Sampietro<sup>1</sup>, L. Dominioni<sup>1</sup>, A. Grande<sup>1</sup>, A. Imperatori<sup>1</sup>, E. d'Amanzio<sup>1</sup>, A. Nicolini<sup>2</sup>  
<sup>1</sup> Center for Thoracic Surgery, University of Insubria, Ospedale di Circolo, Varese, Italy; <sup>2</sup> Respiratory Diseases Unit, Hospital of Sestri Levante, Italy  
5<sup>th</sup> ERCA Lione 2018

**Consecutive use of flow accelerator and cough assistant in patients with neuro-muscular pathology: 2 clinical cases**

E. Mantellini, MD, Pizzorno M.,  
PT Respiratory Rehabilitation Alessandria Hospital Italy  
5<sup>th</sup> ERCA Lione 2018

**Synergic action between expiratory flow accelerator (free aspire) and the in-exsufflator in secretions removal of patient complex: case report**

Giancarlo Garuti MD Respiratory Diseases Dpt;  
Santa Maria Bianca Hospital, Mirandola; Local Health Unit, Modena, Italy  
5<sup>th</sup> ERCA Lione 2018

**Free aspire: a new device for the management of airways clearance in patient with ineffective cough**

L. Bertelli<sup>1</sup>, G. Di Nardo<sup>2</sup>, S. Cazzato<sup>3</sup>, G. Ricci<sup>1</sup>, A. Pession<sup>1</sup>  
<sup>1</sup>Department of Pediatrics, University of Bologna, S. Orsola-Malpighi Hospital; <sup>2</sup>Pediatric Gastroenterology Unit, AORN Santobono-Pausilipon, Naples; <sup>3</sup>Department of Mother and Child Health, Salesi Children's Hospital, Ancona, Italy  
Pediatric Reports 2017; volume 9:7270

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**Evaluation of tolerability, ease of use, perception of mucus encumbrance of an alternative device (FreeAspire ®) for management secretions in patients with bulbar Amyotrophic Lateral Sclerosis**

A. Mattei<sup>1</sup>, E. Rindone<sup>1</sup>, S. Marchisio<sup>1</sup>, M. Bellocchia<sup>1</sup>, R. Ambrosanio<sup>1</sup>, G. Tabbia<sup>1</sup>, E. Clivati<sup>1</sup>, M. Bardessono<sup>1</sup>, C. Bucca<sup>1</sup>  
<sup>1</sup>AOU Città della salute e della Scienza Presidio MOLINETTE SC Pneumologia U - Turin (Italy)  
European Respiratory Journal 2017 50: PA2012; DOI: 10.1183/1393003.congress-2017.PA2012

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**Management of bronchial secretions with Free Aspire in children with cerebral palsy: impact on clinical outcomes and healthcare resources**

G. Garuti<sup>1</sup>, E. Verucchi<sup>2</sup>, I. Fanelli<sup>3</sup>, M. Giovannini<sup>1</sup>, J. C. Winck<sup>4</sup>, M. Lusuardi<sup>5</sup>  
<sup>1</sup>Pneumology Unit, Santa Maria Bianca Hospital, via Fogazzaro 1, Mirandola (MO), 41025 Modena, Italy. <sup>2</sup>Italian Union against Muscular Dystrophy (UILD), Modena and Reggio Emilia Sections, Modena, Italy. <sup>3</sup>public Health Care Department, Child Neuropsychiatric Unit, Modena, Italy. <sup>4</sup>Faculdade de Medicina, Universidade do Porto, Porto, Portugal. <sup>5</sup>Respiratory Rehabilitation Unit, S. Sebastiano Hospital, Correggio (RE), and UILD, Reggio Emilia, Italy.  
Italian Journal of Pediatrics (2016) 42:7

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**Tracheostomized (TCS) patients: Is it possible to manage noninvasively tracheobronchial secretions?**

F. D'Abrosca<sup>1</sup>, I. Prince<sup>1</sup>, S. Belli<sup>1</sup>, C. Carai<sup>1</sup>, M. Moroso<sup>1</sup>, D. Vallese<sup>1</sup>, E. Spada<sup>1</sup>, S. Zaccaria<sup>1</sup>, B. Balbi<sup>1</sup>

<sup>1</sup>Pulmonary Rehabilitation, Fondazione Salvatore Maugeri, I.R.C.C.S., Veruno, Novara, Italy

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**The role of the FREE ASPIRE in the management of the secretions in patients affected by cerebral palsy and its impact in use of health care resources**

G. Garuti<sup>1</sup>, M. Lusuardi<sup>1</sup>, B. Turola<sup>1</sup>, E. Verucchi<sup>2</sup>, M. Giovannini<sup>3</sup>

<sup>1</sup>San Sebastiano Hospital, Correggio (RE), Local Health Authority of Reggio Emilia

<sup>2</sup>Italian Union Against Muscular Dystrophy, Modena Section, Italy

<sup>3</sup>Santa Maria Bianca Hospital, Mirandola (MO), Local Health Authority of Modena

ERCA-JVD Lyon March 2015

**Non-invasive management of tracheobronchial secretions in tracheostomized (TCS) patients**

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ERCA-JVD Lyon March 2015

**Efficacy of Vaküm technology in patients with chronic hypersecretion: a multicentre randomised clinical trial**

E. Venturelli<sup>1,2</sup>, R. Tonelli<sup>2</sup>, D. Righi<sup>1</sup>, V. Furino<sup>1</sup>, P. Cavicchioli<sup>3</sup>, V. Galimberti<sup>4</sup>, A. Albero<sup>5</sup>, A. Colaioli<sup>6</sup>, M. Saleri<sup>7</sup>, C. Castellani<sup>8</sup>, G. Vaghegini<sup>3</sup>, M. Zaurino<sup>4</sup>, P.G. Schiavoni<sup>4</sup>, L. Bianchi<sup>7</sup>, M. Vitacca<sup>7</sup>, L. Palumbo<sup>5</sup>, F. Colombo<sup>5</sup>, F. Pasqua<sup>6</sup>, V. Cardaci<sup>6</sup>, F. Gigliotti<sup>8</sup>, N. Ambrosino<sup>3</sup>, E. Clini<sup>1,2</sup>

<sup>1</sup>Villa Pineta Hospital, Lung Unit and pulmonary Rehabilitation and <sup>2</sup>University of Modena and Reggio Emilia, <sup>3</sup>Rehabilitation Centre Auxilium Vitae- Volterra (Pisa), <sup>4</sup>Pulmonary Rehabilitation Unit of San Giuseppe Multimedica Hospital. Milan, <sup>5</sup>Dpt. Of Rehabilitative Pneumology- Cuasso al Monte (Varese), <sup>6</sup>San Raffaele IRCCS- Rome, <sup>7</sup>Maugeri Foundation IRCCS-Lumezzane (Brescia), <sup>8</sup>Rehabilitation Centredon Gnocchi IRCCS- Florence

ERCA-JVD Lyon March 2015

**FREE ASPIRE versus PEP bottle in the tracheobronchial secretion removal in patients undergoing pulmonary lobectomy: a preliminary prospective and randomized study on 50 patients**

M. D'Andria<sup>1</sup>, E. D'Amazio<sup>1</sup>, L. Cavanna<sup>2</sup>, G. Patrizio<sup>3</sup>, M. Cattoni<sup>4</sup>, C. Sampietro<sup>4</sup>, A. Imperatori<sup>4</sup>, L. Dominioni<sup>4</sup>, N. Rotolo<sup>4</sup>

<sup>1</sup>Degree course in Physiotherapy, School of Medicine, University of Insubria, Varese, Italy; <sup>2</sup>Neuromotor Rehabilitation Unit, Ospedale di Circolo, Varese, Italy; <sup>3</sup>Pneumology Unit, Ospedale di Circolo, Varese, Italy; <sup>4</sup>Center of Thoracic Surgery, University of Insubria, Ospedale di Circolo, Varese, Italy

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#### **Acceptability and effectiveness of a new device in the management of airways clearance**

L. Bertelli MD<sup>1</sup>, S. Cazzato MD<sup>1</sup>, E Di Palma MD<sup>1</sup>, F. Visciotti MD<sup>1</sup>, I. Cecconi MD<sup>2</sup>, A. Pini MD<sup>3</sup>, A. Mazzotta MD<sup>1</sup>, D.M. Cordelli MD<sup>2</sup>, G. Bardasi<sup>1</sup>, A. Pession MD<sup>1</sup>

<sup>1</sup> Pediatrics, University of Bologna, S. Orsola-Malpighi Hospital Bologna, Bologna, Italy; <sup>2</sup> Child Neurology and Psychiatry Unit, S Orsola Malpighi Hospital, University of Bologna, Bologna, Italy; <sup>3</sup> UOC Paediatric Neuropsychiatry, Bellaria-Maggiore Hospital, Bologna, Italy

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#### **FREE-ASPIRE® and chest physiotherapy in bronchial hypersecretion**

P. Cavicchioli, G. Vaghegini, E. Panait, G. Montagnani, D. Berrighi, N. Ambrosino

Pulmonary Rehabilitation Unit - Auxilium Vitae Rehabilitation Center, Volterra, Italy.

UIP 2013 Verona Italy

#### **A non-invasive alternative to suction in tracheostomized (TCS) patients**

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UIP 2013 Verona Italy

#### **One year experience in intensive post-operative rehabilitation (cod.56) in adult after cardiac surgery**

A. Frijia, F. Filippi, A. Gavarini, S. Balloni, M. Ceccarani, I. Luschi, S. Lazzarotti, M. Viscariello, P. Rossi, M. Petrilli  
Fondazione Don Gnocchi, Marina di Massa, Italia

2nd International conference on Respiratory Physiotherapy March 2013 Genova

#### **Free Aspire, an our experience**

A. Bogetti, R. Cara, C. Colliva, Lonati, P. Pisano

2nd International conference on Respiratory Physiotherapy March 2013 Genova

#### **Efficacy of VAKÛM Technology in chronic hypersecretion: a randomised clinical trial**

E. Venturelli<sup>1</sup>, P. Cavicchioli<sup>2</sup>, V. Galimberti<sup>3</sup>, A. Gatta<sup>4</sup>, Antonella Alberò<sup>5</sup>, A. Colaioli<sup>6</sup>, M. Saleri<sup>7</sup>, C. Castellani<sup>8</sup>

<sup>1</sup>Pavullo nel Frignano, MO, Italy; <sup>2</sup>Volterra, PI, Italy; <sup>3</sup>Milano, MI, Italy; <sup>4</sup>Borsalino, AL, Italy; <sup>5</sup>Cuasso al Monte, VA, Italy; <sup>6</sup>Roma, RO, Italy; <sup>7</sup>Lumezzane, BS, Italy; <sup>8</sup>Firenze, FI, Italy

2nd International conference on Respiratory Physiotherapy - Oral presentation - March 2013 Genova

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#### **The effects of the device Free-Aspire in patients with mildly impaired cough effectiveness.**

P. Cavicchioli<sup>1</sup>, G. Vaghegini<sup>1</sup>, M. Ribeiro da Cunha<sup>2</sup>, A. Nicolino<sup>1</sup>

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40° SIMFER – Oral presentation - Sorrento, October 2012

#### **Airways Clearance Techniques in Cystic Fibrosis: Physiology, Devices and the Future**

A. H. Kendrick

Department of Respiratory Medicine, University Hospitals, Bristol

Cystic Fibrosis - Renewed Hopes Through Research - Edited by Dinesh Sriramulu 2012 | cap 22 | 493-518

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**A non-invasive alternative to suction in tracheostomized (TCS) patients.**

F. D'Abrosca, S. Ghizzi, B. Balbi.  
Divisione di Pneumologia Riabilitativa, IRCCS Fondazione Salvatore Maugeri, Veruno, Novara, Italy.  
4<sup>th</sup> ERCA - Barcelona- 2012

**Disobstruction of the proximal airways with Vaküm technology: Study on the efficiency of the Free Aspire system with stage III/IV COPD patients.**

C. Banni Alunno, E. D'amanzio, A. Cabiaglia, L. Palumbo  
Hospital Firm "Ospedale di Circolo" Varese. Italy.  
4<sup>th</sup> ERCA - Barcelona- 2012

**FREE-ASPIRE and chest physiotherapy in bronchial hypersecretion**

P. Cavicchioli, G. Vaghegini, E. Panait, S. Mazzoleni, G. Montagnani, S. Reitano, N. Ambrosino  
Pulmonary Rehabilitation Unit - Auxilium Vitae Rehabilitation Centre, Volterra (PI), The BioRobotics Institute – Scuola Superiore Sant'Anna, Pontedera, Italy.  
4<sup>th</sup> ERCA - Oral presentation -Barcelona- 2012

**Treatment of patients with ineffective cough: observational study about the use of the expiratory flow accelerator FREE ASPIRE**

A. Pedroni <sup>1</sup>, A. Andreoli<sup>2</sup>, A. Cortelezzi<sup>2</sup>  
<sup>1</sup>Graduating in physiotherapy, <sup>2</sup>Luigi Sacco Hospital employees. Milano. Italy.  
4<sup>th</sup> ERCA - Barcelona- 2012

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F. D'Abrosca, B. Balbi  
Il Fisioterapista 5/2010

**The rehabilitative treatment by FREE ASPIRE can reduce pulmonary hypertension in patients with post- cardiac surgery heart failure. Hemodynamic evaluation of echocardiographic parameters**

M. Solinas<sup>1</sup>, G. Mengozzi<sup>2</sup>, A. Di Vincenzo<sup>1</sup>, G. Solarino<sup>2</sup>, A. Frijia<sup>1</sup>, M- Petrilli<sup>1</sup>.  
<sup>1</sup>O.U. of Cardiovascular and Respiratory Rehabilitation, Don Gnocchi-ONLUS Foundation, Marina di Massa – Italy; <sup>2</sup> Cardiovascular Department Hospital University of Pisa  
ERS 2010 Barcellona

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Medical Products Research Srl