## Lavori scientifici / Clinical evidences

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### **Comparison of two mechanical insufflation-exsufflation devices in patients with amyotrophic lateral sclerosis: a preliminary study**

Antonello Nicolini<sup>1</sup>, Paola Prato<sup>2</sup>, Laura Beccarelli<sup>2</sup>, Bruna Grecchi<sup>3</sup>, Giancarlo Garuti<sup>4</sup>, Paolo Banfi<sup>5</sup>, Francesco D'Abrosca<sup>6</sup>

<sup>1</sup>Respiratory Diseases Unit, General Hospital, Sestri Levante, Genova, Italy - antonellonicolini@gmail.com.

<sup>2</sup>Respiratory Diseases Unit, General Hospital, Sestri Levante, Genova, Italy.

<sup>3</sup>Physical Medicine and Rehabilitation, General Hospital, Sestri Levante, Genova, Italy.

<sup>4</sup>Department of Pulmonology, General Hospital, Mirandola, Modena, Italy.

<sup>5</sup>Pulmonary Rehabilitation Unit, Don Gnocchi Foundation, IRCCS, Milan, Italy.

<sup>6</sup>ARIR, Italian Association of Respiratory Physiotherapists.

Panminerva Med. 2021 Mar 9. doi: 10.23736/S0031-0808.21.04104-5. Epub ahead of print. PMID: 33687180.

### **EFA® technology efficacy on oropharyngeal secretions management in subjects with DOC**

Alessandra Cima<sup>1</sup>, Manuela Locatelli<sup>2</sup>, Francesco D'Abrosca<sup>3,4</sup>

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<sup>2</sup>ICS Maugeri – Lissone Italy

<sup>3</sup>Italian Association of Respiratory Physiotherapists (ARIR), Milan, Italy

<sup>4</sup>Department of Health Sciences - University of Eastern Piedmont

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### **Comparison of two mechanical insufflation-exsufflation devices in patients with amyotrophic lateral sclerosis: a preliminary study**

Nicolini Antonello<sup>1</sup>, Prato Paola<sup>1</sup>, Beccarelli Laura<sup>1</sup>, Grecchi Bruna<sup>2</sup>, Garuti Giancarlo<sup>3</sup>, Banfi Paolo<sup>4</sup>, Francesco D'Abrosca<sup>5</sup>

<sup>1</sup>Respiratory Diseases Unit, General Hospital, Sestri Levante, Italy

<sup>2</sup>Physical Medicine and Rehabilitation, General Hospital, Sestri Levante, Italy 3.Department of Pulmonology, General Hospital, Mirandola, Italy

<sup>4</sup>Pulmonary rehabilitation Unit, Don Gnocchi foundation, IRCCS, Milan

<sup>5</sup>ARIR, Italian Association of Respiratory Physiotherapists.

Panminerva Med. 2021 Mar 9. doi: 10.23736/S0031-0808.21.04104-5. Epub ahead of print. PMID: 33687180.

### **Comparison of an expiratory flow accelerator device versus positive expiratory pressure for tracheobronchial airway clearance after lung cancer lobectomy: a preliminary study**

Nicola Rotolo<sup>1</sup>, Maria Cattoni<sup>1</sup>, Michele D'Andria<sup>1</sup>, Laura Cavanna<sup>1</sup>, Giorgia Patrizio<sup>1</sup>, Andrea Imperatori<sup>1</sup>, Antonello Nicolini<sup>2</sup>

<sup>1</sup>Center for Thoracic Surgery, University of Insubria, Ospedale di Circolo, Varese, Italy

<sup>2</sup>Respiratory Diseases Unit, Hospital of Sestri Levante, Italy

Physiotherapy. 2021 Mar;110:34-41. doi: 10.1016/j.physio.2019.01.011. Epub 2019 Jan 26. PMID: 33563372.

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**Effectiveness of expiratory flow acceleration in patients with Parkinson's disease and swallowing deficiency: A preliminary study**

Giulio Riboldazzi<sup>1</sup>, Giada Spinazza<sup>2</sup>, Laura Beccarelli<sup>3</sup>, Paola Prato<sup>3</sup>, Bruna Grecchi<sup>4</sup>, Francesco D'Abrosca<sup>5</sup>, Antonello Nicolini<sup>3</sup>

<sup>1</sup> Parkinson's Disease Unit, Gaetano and Piera Borghi Foundation, Brebbia, Italy

<sup>2</sup> Neurology Rehabilitation Department, Hospital of Cuasso al Monte, Italy

<sup>3</sup> Pulmonary Rehabilitation Unit, General Hospital, Sestri Levante, Italy

<sup>4</sup> Physical Medicine and Rehabilitation Department, General Hospital, Sestri Levante, Italy

<sup>5</sup> Italian Association of Respiratory Physiotherapists (ARIR), Milan, Italy

Clinical Neurology and Neurosurgery 199 (2020) 106249

**Expiratory Flow Accelerator (EFA®): home-based airway clearance daily program in systemic sclerosis patients with mucus hypersecretion, productive cough and dyspnoea. A pilot study.**

Silvia Faverzani<sup>1</sup>, A. Ariani<sup>1</sup>, A. Becciolini<sup>1</sup>, E. Di Donato<sup>1</sup>, F. Nocera<sup>1</sup>, F. Mozzani<sup>1</sup>, D. Santilli<sup>1</sup>, M. Riva<sup>1</sup>, L. Monica<sup>1</sup>, A. Barbieri<sup>1</sup>, L. Barone<sup>1</sup>, V. Alfieri<sup>1</sup>, A. Frizzelli<sup>1</sup>, A. Chetta<sup>1</sup>, E. Crisafulli<sup>2</sup>, M. Marvisi<sup>3</sup>

<sup>1</sup>University Hospital of Parma, Italy

<sup>2</sup>University Hospital of Verona, Italy

<sup>3</sup>Figlie di S. Camillo Cremona, Italy

EULAR 2020

**Airway Clearance with Expiratory Flow Accelerator Technology: Effectiveness of the "Free Aspire" Device in Patients with Severe COPD.**

Giorgia Patrizio<sup>1</sup>, Michele D'Andria<sup>2</sup>, Francesco D'Abrosca<sup>3</sup>, Antonella Cabiaglia<sup>1</sup>, Fabio Tanzi<sup>4</sup>, Giancarlo Garuti<sup>5</sup>, Antonello Nicolini<sup>6</sup>

<sup>1</sup>Pulmonology and Rehabilitation Unit, Cuasso al Monte Hospital, Varese, Italy

<sup>2</sup>Physiotherapy Degree Course, Insubria University, Varese, Italy

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<sup>4</sup>Physical Medicine and Rehabilitation, General Hospital, Varese, Italy

<sup>5</sup>Pulmonology Unit, Santa Maria Bianca Hospital, Mirandola, Italy

<sup>6</sup>Respiratory Diseases Unit, General Hospital Sestri Levante, Sestri Levante, Italy

Turkish Thoracic Journal 2019; 20(4): 209-15

**Expiratory Flow Accelerator (EFA) technique on mucus hypersecretion of COPD patients with reduced cough efficiency after a severe exacerbation.**

Elisabetta Zampogna<sup>1</sup>, Ernesto Crisafulli<sup>2,3</sup>, Michele D'Andria<sup>4</sup>, Cristina Gregorini<sup>4</sup>, Giorgio Bellelli<sup>1</sup>, Etienne Lucini<sup>1</sup>, Silvia Faverzani<sup>5</sup>, Dina Visca<sup>6</sup>, Antonio Spanevello<sup>6</sup>, Nicola Schiavone<sup>7</sup>, Alfredo Chetta<sup>5</sup>, Andrea Zanini<sup>7</sup>

<sup>1</sup>Division of Pulmonary Rehabilitation, Istituti Clinici Scientifici Maugeri, IRCCS, Tradate, Italy

<sup>2</sup>Department of Medicine, Respiratory Medicine Unit, University of Verona and Azienda Ospedaliera Universitaria Integrata of Verona, Verona, Italy

<sup>3</sup>Department of Medicine, Section of Internal Medicine, University of Verona and Azienda Ospedaliera Universitaria Integrata of Verona, Verona, Italy

<sup>4</sup>Division of General Medicine, Ospedale Malcantonese, Castelrotto, Switzerland

<sup>5</sup>Department of Medicine and Surgery, Respiratory Disease and Lung Function Unit, University of Parma, Parma, Italy

<sup>6</sup>Department of Medicine and Surgery, Respiratory Diseases, University of Insubria, Varese-Como, Italy

<sup>7</sup>Pulmonary Rehabilitation, Clinic of Rehabilitation, Ente Ospedaliero Cantonale, Novaggio, Switzerland

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### **A pilot study on the non-invasive management of tracheobronchial secretions in tracheostomised patients**

S. Belli<sup>1</sup>, D. Cattaneo<sup>1</sup>, F. D'Abrosca<sup>2</sup>, I. Prince<sup>1</sup>, G. Savio<sup>1</sup>, B. Balbi<sup>1</sup>

<sup>1</sup>Istituti Clinici Scientifici Maugeri IRCCS, Respiratory Rehabilitation Unit of Veruno Institute, Veruno, Italy

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Clinical Respiratory Journal 2019; 13.10: 637-642

### **Airway clearance management with vakum technology in subjects with ineffective cough: A pilot study on the efficacy, acceptability evaluation, and perception in children with cerebral palsy.**

L. Bertelli<sup>1</sup>, G. Bardasi<sup>1</sup>, S. Cazzato<sup>2</sup>, E. Di Palma<sup>1</sup>, M. Gallucci<sup>1</sup>, G. Ricci<sup>1</sup>, A. Pession<sup>1</sup>

<sup>1</sup>Department of Pediatrics, University of Bologna, S. Orsola-Malpighi Hospital, Bologna, Italy; <sup>2</sup>Department of Mother and Child Health, Salesi Children's Hospital, Ancona Italy.

Pediatric Allergy, Immunology, and Pulmonology, 2019, 32.1: 23-27

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### **KALOS device in SMA patients: a case report**

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### **Management of the airway in the tracheostomy patient in critical areas: a comparison of two devices. A pilot study.**

C. Blardone<sup>1</sup>, L. Cavanna<sup>2</sup>, M. Bertoni<sup>2</sup>, P. Severgnini<sup>3</sup>

<sup>1</sup>Pulmonology Department ASST-Sette Laghi, Varese, Italy; <sup>2</sup>Rehabilitation Department ASST-Sette Laghi, Varese, Italy; <sup>3</sup>University of Insubria -DBSV -ICU ASST Sette Laghi Varese, Italy

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### **Short-term effectiveness of the expiratory flow accelerator (EFA) on the polmonary hyperinsufflation: a case report.**

S. Musella, E. Sciarillo, F. Stefanelli, A. Di Giorgio, M. Lanza, R. Cauteruccio, A. Annunziata, I. Meoli

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### **Validation study of a protocol for airway clearance in patients with ineffective cough**

S. Belli, D. Cattaneo, F. Masocco, M. Peroni, G. Maggioni, C. Marcassa, B. Bruno

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### **Effectiveness of an expiratory flow accelerator device (Free Aspire) in dysphagic patients with Parkinson's disease: a preliminary study**

A. Nicolini, B. Grecchi

Respiratory Diseases Unit, General Hospital, Sestri Levante, Italy

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63° Congresso Nazionale SIGG – Roma 2018

**The kalos device: experimental protocol for airway clearance in patients with ineffective cough**

S. Belli<sup>1</sup>, F. Masocco<sup>1</sup>, D. Cattaneo<sup>1</sup>, G. Savio<sup>1</sup>, B. Balbi<sup>1</sup>

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5<sup>th</sup> ERCA Lione 2018

**A pilot study on the non-invasive management of tracheobronchial secretions in tracheotomised patients**

S. Belli<sup>1</sup>, D. Cattaneo<sup>1</sup>, F. D'Abrosca<sup>2</sup>, I. Prince<sup>1</sup>, G. Savio<sup>1</sup>, B. Balbi<sup>1</sup>.

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5<sup>th</sup> ERCA Lione 2018

**Management of the airway in the tracheostomy patient in critical areas: a comparison of two devices for a pilot study.**

L. Cavanna<sup>1</sup> PT, C. Blardone<sup>2</sup> PT, M. Bertoni MD<sup>3</sup>, P. Severgnini MD<sup>4</sup>

<sup>1,3</sup> Rehabilitation Department ASST-Sette Laghi, Varese; <sup>2</sup> Pulmonology Department ASST-Sette Laghi, Varese; <sup>4</sup> University of Insubria - DBSV – ICU ASST Settelaghi Varese, Italy;

5<sup>th</sup> ERCA Lione 2018

**Running title: EFA and pep physiotherapy after pulmonary lobectomy**

N. Rotolo<sup>1</sup>, M. Cattoni<sup>1</sup>, M. D'Andria<sup>1</sup>, L. Cavanna<sup>1</sup>, G. Patrizio<sup>1</sup>, C. Sampietro<sup>1</sup>, L. Dominioni<sup>1</sup>, A. Grande<sup>1</sup>, A. Imperatori<sup>1</sup>, E. d'Amanzio<sup>1</sup>, A. Nicolini<sup>2</sup>

<sup>1</sup> Center for Thoracic Surgery, University of Insubria, Ospedale di Circolo, Varese, Italy; <sup>2</sup> Respiratory Diseases Unit, Hospital of Sestri Levante, Italy

5<sup>th</sup> ERCA Lione 2018

**Consecutive use of flow accelerator and cough assistant in patients with neuro-muscular pathology: 2 clinical cases**

E. Mantellini, MD, Pizzorno M.,

PT Respiratory Rehabilitation Alessandria Hospital Italy

5<sup>th</sup> ERCA Lione 2018

**Synergic action between expiratory flow accelerator (free aspire) and the in-exsufflator in secretions removal of patient complex: case report**

Giancarlo Garuti MD Respiratory Diseases Dpt;

Santa Maria Bianca Hospital, Mirandola; Local Health Unit, Modena, Italy

5<sup>th</sup> ERCA Lione 2018

**Free aspire: a new device for the management of airways clearance in patient with ineffective cough**

L. Bertelli<sup>1</sup>, G. Di Nardo<sup>2</sup>, S. Cazzato<sup>3</sup>, G. Ricci<sup>1</sup>, A. Pession<sup>1</sup>

<sup>1</sup>Department of Pediatrics, University of Bologna, S. Orsola-Malpighi Hospital; <sup>2</sup>Pediatric Gastroenterology Unit, AORN Santobono-Pausilipon, Naples; <sup>3</sup>Department of Mother and Child Health, Salesi Children's Hospital, Ancona, Italy

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XXXVI SIAF 2017 Siena

**Evaluation of tolerability, ease of use, perception of mucus encumbrance of an alternative device (FreeAspire<sup>®</sup>) for management secretions in patients with bulbar Amyotrophic Lateral Sclerosis**

A. Mattei<sup>1</sup>, E. Rindone<sup>1</sup>, S. Marchisio<sup>1</sup>, M. Bellocchia<sup>1</sup>, R. Ambrosanio<sup>1</sup>, G. Tabbia<sup>1</sup>, E. Clivati<sup>1</sup>, M. Bardesson<sup>1</sup>, C. Bucca<sup>1</sup>

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**Management of bronchial secretions with Free Aspire in children with cerebral palsy: impact on clinical outcomes and healthcare resources**

G. Garuti<sup>1</sup>, E. Verucchi<sup>2</sup>, I. Fanelli<sup>3</sup>, M. Giovannini<sup>1</sup>, J. C. Winck<sup>4</sup>, M. Lusuardi<sup>5</sup>

<sup>1</sup>Pneumology Unit, Santa Maria Bianca Hospital, via Fogazzaro 1, Mirandola (MO), 41025 Modena, Italy. <sup>2</sup>Italian Union against Muscular Dystrophy (UILDM), Modena and Reggio Emilia Sections, Modena, Italy. <sup>3</sup>public Health Care Department, Child Neuropsychiatric Unit, Modena, Italy. <sup>4</sup>Faculdade de Medicina, Universidade do Porto, Porto, Portugal. <sup>5</sup>Respiratory Rehabilitation Unit, S. Sebastiano Hospital, Correggio (RE), and UILDM, Reggio Emilia, Italy.

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**Tracheostomized (TCS) patients: Is it possible to manage noninvasively tracheobronchial secretions?**

F. D'Abrosca<sup>1</sup>, I. Prince<sup>1</sup>, S. Belli<sup>1</sup>, C. Carai<sup>1</sup>, M. Moroso<sup>1</sup>, D. Vallese<sup>1</sup>, E. Spada<sup>1</sup>, S. Zaccaria<sup>1</sup>, B. Balbi<sup>1</sup>

<sup>1</sup>Pulmonary Rehabilitation, Fondazione Salvatore Maugeri, I.R.C.C.S., Veruno, Novara, Italy

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**The role of the FREE ASPIRE in the management of the secretions in patients affected by cerebral palsy and its impact in use of health care resources**

G. Garuti<sup>1</sup>, M. Lusuardi<sup>1</sup>, B. Turola<sup>1</sup>, E. Verucchi<sup>2</sup>, M. Giovannini<sup>3</sup>

<sup>1</sup>San Sebastiano Hospital, Correggio (RE), Local Health Authority of Reggio Emilia

<sup>2</sup>Italian Union Against Muscular Dystrophy, Modena Section, Italy

<sup>3</sup>Santa Maria Bianca Hospital, Mirandola (MO), Local Health Authority of Modena

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### **Non-invasive management of tracheobronchial secretions in tracheostomized (TCS) patients**

F. D'Abrosca, S. Belli, C. Carai, M. Moroso, I. Pince, D. Vallese, B. Balbi  
*Divisione di Pneumologia Riabilitativa, Fondazione S. Maugeri, IRCCS, Veruno (NO), Italy*  
ERCA-JVD Lyon March 2015

### **Efficacy of Vaküm technology in patients with chronic hypersecretion: a multicentre randomised clinical trial**

E. Venturelli<sup>1,2</sup>, R. Tonelli<sup>2</sup>, D. Righi<sup>1</sup>, V. Furino<sup>1</sup>, P. Cavicchioli<sup>3</sup>, V. Galimberti<sup>4</sup>, A. Albero<sup>5</sup>, A. Colaioli<sup>6</sup>, M. Saleri<sup>7</sup>, C. Castellani<sup>8</sup>, G. Vaghegginì<sup>3</sup>, M. Zaurino<sup>4</sup>, P.G. Schiavoni<sup>4</sup>, L. Bianchi<sup>7</sup>, M. Vitacca<sup>7</sup>, L. Palumbo<sup>5</sup>, F. Colombo<sup>5</sup>, F. Pasqua<sup>6</sup>, V. Cardaci<sup>6</sup>, F. Gigliotti<sup>8</sup>, N. Ambrosino<sup>3</sup>, E. Clini<sup>1,2</sup>

<sup>1</sup>Villa Pineta Hospital, Lung Unit and pulmonary Rehabilitation and <sup>2</sup>University of Modena and Reggio Emilia, <sup>3</sup>Rehabilitation Centre Auxilium Vitae- Volterra (Pisa), <sup>4</sup>Pulmonary Rehabilitation Unit of San Giuseppe Multimedica Hospital. Milan, <sup>5</sup>Dpt. Of Rehabilitative Pneumology- Cuasso al Monte (Varese), <sup>6</sup>San Raffaele IRCCS- Rome, <sup>7</sup>Maugeri Foundation IRCCS-Lumezzane (Brescia), <sup>8</sup>Rehabilitation Centredon Gnocchi IRCCS- Florence

ERCA-JVD Lyon March 2015

### **FREE ASPIRE versus PEP bottle in the tracheobronchial secretion removal in patients undergoing pulmonary lobectomy: a preliminary prospective and randomized study on 50 patients**

M. D'Andria<sup>1</sup>, E. D'Amancio<sup>1</sup>, L. Cavanna<sup>2</sup>, G. Patrizio<sup>3</sup>, M. Cattoni<sup>4</sup>, C. Sampietro<sup>4</sup>, A. Imperatori<sup>4</sup>, L. Dominioni<sup>4</sup>, N. Rotolo<sup>4</sup>  
<sup>1</sup>Degree course in Physiotherapy, School of Medicine, University of Insubria, Varese, Italy; <sup>2</sup>Neuromotor Rehabilitation Unit, Ospedale di Circolo, Varese, Italy; <sup>3</sup>Pneumology Unit, Ospedale di Circolo, Varese, Italy; <sup>4</sup>Center of Thoracic Surgery, University of Insubria, Ospedale di Circolo, Varese, Italy

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### **Acceptability and effectiveness of a new device in the management of airways clearance**

L. Bertelli MD<sup>1</sup>, S. Cazzato MD<sup>1</sup>, E Di Palmo MD<sup>1</sup>, F. Visciotti MD<sup>1</sup>, I. Cecconi MD<sup>2</sup>, A. Pini MD<sup>3</sup>, A. Mazzotta MD<sup>1</sup>, D.M. Cordelli MD<sup>2</sup>, G. Bardasi<sup>1</sup>, A. Pession MD<sup>1</sup>

<sup>1</sup> Pediatrics, University of Bologna, S. Orsola-Malpighi Hospital Bologna, Bologna, Italy; <sup>2</sup> Child Neurology and Psychiatry Unit, S Orsola Malpighi Hospital, University of Bologna, Bologna, Italy; <sup>3</sup> UOC Paediatric Neuropsychiatry, Bellaria-Maggiore Hospital, Bologna, Italy

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### **FREE-ASPIRE® and chest physiotherapy in bronchial hypersecretion**

P. Cavicchioli, G. Vaghegginì, E. Panait, G. Montagnani, D. Berrighi, N. Ambrosino  
*Pulmonary Rehabilitation Unit - Auxilium Vitae Rehabilitation Center, Volterra, Italy.*  
UIP 2013 Verona Italy

Medical Products Research Srl

Sede legale: P.zza San Magno, 7 - 20025 Legnano (MI)  
Stabilimento: Via Novara, 68 - 20025 Legnano (MI)  
Tel. +39.0331.597992 Fax +39.0331.485089  
<http://www.mpr-italy.it> e-mail: [info@mpr-italy.it](mailto:info@mpr-italy.it)

Capitale sociale €100.000 int.ver. - C.C.I.A.A. Milano n° 04984420960 - Rea Milano n° 1787932  
Codice fiscale 04984420960 Partita IVA 04984420960

#### **A non-invasive alternative to suction in tracheostomized (TCS) patients**

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UIP 2013 Verona Italy

#### **One year experience in intensive post-operative rehabilitation (cod.56) in adult after cardiac surgery**

A. Frijia, F. Filippi, A. Gavarini, S. Balloni, M. Ceccarani, I. Luschi, S. Lazzarotti, M. Viscariello, P. Rossi, M. Petrilli

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2nd International conference on Respiratory Physiotherapy March 2013 Genova

#### **Free Aspire, an our experience**

A. Bogetti, R. Cara, C. Colliva, Lonati, P. Pisano

2nd International conference on Respiratory Physiotherapy March 2013 Genova

#### **Efficacy of VAKÜM Technology in chronic hypersecretion: a randomised clinical trial**

E. Venturelli<sup>1</sup>, P. Cavicchioli<sup>2</sup>, V. Galimberti<sup>3</sup>, A. Gatta<sup>4</sup>, Antonella Albero<sup>5</sup>, A. Colaioli<sup>6</sup>, M. Saleri<sup>7</sup>, C. Castellani<sup>8</sup>

<sup>1</sup>Pavullo nel Frignano, MO, Italy; <sup>2</sup>Volterra, PI, Italy; <sup>3</sup>Milano, MI, Italy; <sup>4</sup>Borsalino, AL, Italy; <sup>5</sup>Cuasso al Monte, VA, Italy; <sup>6</sup>Roma, RO, Italy; <sup>7</sup>Lumezzane, BS, Italy; <sup>8</sup>Firenze, FI, Italy

2nd International conference on Respiratory Physiotherapy - Oral presentation - March 2013 Genova

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#### **The effects of the device Free-Aspire in patients with mildly impaired cough effectiveness.**

P. Cavicchioli<sup>1</sup>, G. Vaghegini<sup>1</sup>, M. Ribeiro da Cunha<sup>2</sup>, A. Nicolino<sup>1</sup>

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40° SIMFER – Oral presentation - Sorrento, October 2012

#### **Airways Clearance Techniques in Cystic Fibrosis: Physiology, Devices and the Future**

A. H. Kendrick

*Department of Respiratory Medicine, University Hospitals, Bristol*

Cystic Fibrosis - Renewed Hopes Through Research - Edited by Dinesh Sriramulu 2012 |cap 22| 493-518

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SIMRI 2012 Genova

#### **A non-invasive alternative to suction in tracheostomized (TCS) patients.**

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4<sup>th</sup> ERCA - Barcelona- 2012

#### **Disobstruction of the proximal airways with Vaküm technology: Study on the efficiency of the Free Aspire system with stage III/IV COPD patients.**

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4<sup>th</sup> ERCA - Barcelona- 2012

#### **FREE-ASPIRE and chest physiotherapy in bronchial hypersecretion**

P. Cavicchioli, G. Vaghegini, E. Panait, S. Mazzoleni, G. Montagnani, S. Reitano, N. Ambrosino

*Pulmonary Rehabilitation Unit - Auxilium Vitae Rehabilitation Centre, Volterra (PI), The BioRobotics Institute – Scuola Superiore Sant'Anna, Pontedera, Italy.*

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**Treatment of patients with ineffective cough: observational study about the use of the expiratory flow accelerator FREE ASPIRE**

A. Pedroni<sup>1</sup>, A. Andreoli<sup>2</sup>, A. Cortelezzi<sup>2</sup>

<sup>1</sup>Graduating in physiotherapy, <sup>2</sup>Luigi Sacco Hospital employees. Milano. Italy.

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**The rehabilitative treatment by FREE ASPIRE can reduce pulmonary hypertension in patients with post- cardiac surgery heart failure. Hemodynamic evaluation of echocardiographic parameters**

M. Solinas<sup>1</sup>, G. Mengozzi<sup>2</sup>, A. Di Vincenzo<sup>1</sup>, G. Solarino<sup>2</sup>, A. Frijia<sup>1</sup>, M- Petrilli<sup>1</sup>.

<sup>1</sup>O.U. of Cardiovascular and Respiratory Rehabilitation, Don Gnocchi-ONLUS Foundation, Marina di Massa – Italy; <sup>2</sup> Cardiovascular Department Hospital University of Pisa

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